

How to measure a bariatric patient

Shapes can be difficult to measure. Measure the person - not the wheelchair.

	11 35
Seat Depth Measure seat depth in both sides. Can the patient use legs for active movements? Then make sure to calculate it in the seat depth. There has to be enough space, so the patient can flex/extend the knee joint to create active "walking" movements in the wheelchair. Inches/cm:	
Full back height Measure full back height. Inches/cm:	
Lower back support Measure lower back support. Measure from the edge of the seat to the inferior angle of scapula. Inches/cm:	
Gluteal shelf height Measure from the edge of the seat to the top of the gluteal shelf. Inches/cm:	



Gluteal shelf depth The depth may vary if there is a pelvic rotation. Inches/cm:	
Shoulder width Make sure the patient sits as straight as possible. Inches/cm:	
Abdominal width Measure abdominal width but be aware of abdominal asymmetry. Inches/cm: Is the abdominal tissue hard or soft? Hard Soft Hips width Inches/cm:	
Distance between knees Measure on the lateral side of the knees. Inches/cm:	



Lower leg length Does the patient wear shoes or not? Inches/cm with/without shoes:	
Ankle joint Can the ankle joint reach neutral position? Inversion Eversion Increased plantar flexion Increased dorsal flexion	
90° over the hip joint Can the patient, passive or active, reach 90° over the hip joint? Note the range of motion of the hip joint. Right, degrees: Left, degrees:	90°
Knee flexion Note the range of motion. Right, degrees: Left, degrees:	90°
Headrest Does the patient have full head and neck control? Yes No If not what kind of headrest does the patient need? Center headrest Lateral headrest	





Patient transfer How is patient transfers performed? Self-independent (no equipment) With standing transfer products Sideways over armrest Hoist	
Propelling the wheelchair How does the patient propel the wheelchair? With arms With legs With arms and legs Need of electric propulsion	
Armrest support While seated and during transfers? Full support Elbow support	
Weight In kilos or lbs:	

